american diagnostica inc.

500 West Avenue, Stamford, CT 06902
Tel. (203) 602-7777 • Fax. (203) 602-2221

ENGLISH

IMUBIND[®] ADAMTS13 Autoantibody ELISA

REF 814

american diagnostica GmbH
Kaplaneigasse 35, D-64319 Pfungstadt, Germany

INTENDED USE

The IMUBIND[®] ADAMTS13 Autobody ELISA is intended for the measurement of ADAMTS13 IgG autoantibodies in human plasma. The assay is limited to "For Research Use Only" in the United States.

EXPLANATION OF THE TEST

ADAMTS13, also known as von Willebrand Factor (vWF) cleaving protease, is a zinc metalloproteinase that cleaves ultra large vWF multimers (UL-vWF) at the Tyr(1605)-Met(1606) bond located in the A2 region of vWF.¹ Studies have shown that low levels of ADAMTS13 activity are associated with Thrombotic Thrombocytopenia Purpura (TTP), a life-threatening hematological condition characterized by low platelet count, microvascular thrombi, red cell fragmentation, CNS and renal complications.^{2,3} A deficiency or low level of ADAMTS13 activity (<5%) may lead to an accumulation of UL-vWF multimers.⁴ The UL-vWF multimers will bind to receptors on platelets inducing platelet aggregation and formation of intravascular thrombi.

Congenital TTP is a rare inheritable disorder resulting from mutations in the ADAMTS13 gene. Genetic alterations have been identified at many sites within the ADAMTS13 gene and result in the production of non-functional ADAMTS13 protein.^{5,6} The acquired form of TTP is an autoimmune-like disorder caused by the development of autoantibodies to ADAMTS13 that inhibit its activity.⁷ The development of ADAMTS13 autoantibodies has been associated with certain drugs such as ticlopidine, clopidogrel, quinine and cyclosporine. The measurement of ADAMTS13 autoantibodies may provide information toward the diagnosis and management of patients with acquired TTP.

PRINCIPLE OF THE METHOD

Diluted plasma samples are added to microwells coated with a full-length recombinant ADAMTS13. During an incubation period, autoantibodies to ADAMTS13 present in the sample will bind to the protein coated to the wells. Following a washing step, a goat anti-human IgG antibody labeled with horseradish peroxidase (HRP) is added to the microwells and binds to IgG type ADAMTS13 autoantibodies that have bound to the coated microwells. Following another washing step, the addition of a perborate-3,3',5,5'-tetramethylbenzidine (TMB) substrate and its subsequent reaction with the HRP present generates a blue colored solution. The reaction is stopped by adding sulfuric acid, which turns the solution yellow. Measuring the solution absorbance at 450 nm and comparing the value with those of a standard curve determines the level of ADAMTS13 autoantibodies in the diluted plasma sample. ADAMTS13 autoantibody levels are expressed in Arbitrary Units/mL (AU/mL). One AU/mL is equal to 1 µg/mL of affinity purified human anti-ADAMTS13 IgG.

REAGENTS

- R1 ADAMTS13 Coated Microwells, 96, with acetate cover sheet
- R2 Assay Buffer, 15 mL, 2 vials (lyophilized)
- R3 Plasma Standard, 60 AU/mL, 0.45mL, 2 vials (lyophilized)
- R4 Positive Control, 0.35 mL, 2 vials (lyophilized)
- R5 Detection Antibody, goat anti-human IgG-HRP, 140 µL
- R6 Substrate, TMB, 11 mL
- R7 Wash Buffer, PBS w/0.05% Tween 20, pH 7.4, 1 Liter (lyophilized)

WARNING

Source material for some of the reagents in this kit is of human origin. This material has been found to be non-reactive for Hepatitis B Surface Antigen (HBsAg), Hepatitis C Virus (HCV) and Human Immunodeficiency Virus Type 1 and Type 2 (HIV-1, HIV-2) using FDA approved methods. As no known test method provides complete assurance that products derived from human blood will not transmit HBsAg, HCV, HIV-1, HIV-2 or other blood-borne pathogens, reagents should be handled as recommended for any potentially infectious human specimen. Discard all waste associated with test specimens and human source reagents in a biohazard waste container.

Limited for research use only in the United States. For *in vitro* use only. Not for internal use in humans or animals. Do not use the kit components beyond the stated expiration date. Do not mix reagents from different kits. Avoid microbial contamination of the reagents. Do not smoke, eat or drink in areas in which specimens or kit reagents are handled. Do not pipette reagents by mouth. Wear laboratory coat and disposable gloves throughout the test procedure and wash hands thoroughly afterwards. Avoid splashing or aerosol formation.

REAGENT PREPARATION AND STORAGE

Unopened and lyophilized reagents are stable until the expiration date printed on the box when stored as instructed.

1. **96 ADAMTS13 Coated Microwells, 6 x 16 well strips:** Once removed from the foil pouch, the microwells must be used within 30 minutes. Unused microwells may be stored at 2°-8°C for 4 weeks when sealed in the original pouch with the desiccant present, protected from any moisture.
2. **Assay Buffer:** Reconstitute with 15 mL of filtered deionized/distilled water. Gently mix the contents of the vial. Assay Buffer may be used for up to 4 weeks when stored at 2°-8°C.
3. **Plasma Standard:** Reconstitute with 0.45 mL of Assay Buffer. Standard may be aliquoted and stored at -20°C for up to 6 months.

4. **Positive Control:** Reconstitute with 0.35 mL of Assay Buffer. Once reconstituted, the control may be stored for up to 4 weeks at -20°C. The control must be diluted 1:2 or 1:4 with Assay Buffer prior to assaying.
5. **Detection Antibody:** Prepare working strength Detection Antibody by diluting 1:100 with Assay Buffer. For using all 96 microwells at one time, dilute 100 µL of Detection Antibody Conjugate up to 10 mL in Assay Buffer. For using less than 96 microwells, dilute 20 µL of Detection Antibody up to 2 mL in Assay Buffer for each 16 microwell strip that will be used. Working strength Detection Antibody may be stored in the dark for 4 hours at room temperature or for 72 hours at 2°-8°C.
6. **Substrate, TMB:** Supplied ready to use. It may be used until the expiration date stated on the vial when stored in the dark at 2°-8°C.
7. **Wash Buffer:** Dissolve the contents of the packet with 1 Liter of filtered deionized/distilled water. After reconstitution, the Wash Buffer may be used for up to 4 weeks when stored at 2°-8°C.

SPECIMEN COLLECTION AND PREPARATION

Only citrate collected platelet poor plasma may be used for this assay. Do Not Use EDTA. See "Collection, Transport and Processing of Blood Specimens for Testing Plasma-based Coagulation Assays; Approved Guidelines-Fourth Edition", NCCLS Document H21-A4, Vol. 23, No. 35, December 2003. Plasma collection should be performed as follows:

1. Collect 9 parts of blood into 1 part of 3.2% (0.109 M) trisodium citrate anticoagulant solution.
2. Centrifuge the blood sample at 10,000 x g for 15 minutes.
3. Plasma should be stored at 2°-8°C and assayed within 4 hours. Alternatively, plasma may be stored at -20°C for up to 6 months.
4. Frozen plasma should be thawed rapidly at 37°C. Thawed plasmas should be stored at 2°-8°C and assayed within 4 hours.

PROCEDURE

MATERIALS PROVIDED – SEE REAGENTS

MATERIAL REQUIRED BUT NOT PROVIDED

0.22 µm filtered deionized H₂O
50-300 µL eight channel multi-pipette
0-200 µL, 200-1000 µL single pipettes
microwell plate reader for reading absorbance at 450 nm
microwell plate orbital shaker, microwell plate washer (optional)
0.5 M H₂SO₄ **Caution:** Handle Sulphuric acid with great care. Avoid any skin and eye contact. Wear protection glasses and gloves when handling.

Standard Preparation

1. Open the foil pouch and remove the microwell strips/frame assembly. Remove the strips that will not be used, return them to the foil pouch and tightly reseal the pouch with the desiccant inside. Store the foil pouch at 2 - 8°C.
2. Reconstitute the Plasma Standard as instructed under REAGENT PREPARATION. Prepare five (5) serial dilutions of the standard as follows. Add 200 µL of Standard to microwells A1/A2. Add 100 µL of Assay Buffer to microwells B1/B2 – F1/F2.
3. Serially dilute the standard by pipetting 100 µL of the standard from microwells A1/A2 into microwells B1/B2. Mix and pipette 100 µL from wells B1/B2 to wells C1/C2. Repeat this process through wells F1/F2. Remove and discard 100 µL from wells F1/F2. The concentrations of the serially diluted standard will be 60, 30, 15, 7.5, 3.75 and 1.87 AU/mL respectively. Add 100 µL of Assay Buffer to wells G1/G2 to serve as the 0 AU/mL.

Assay Procedure

4. Dilute each plasma sample 1:20 (1 part plasma + 19 parts Assay Buffer). Add 100 µL of Positive Control or diluted sample to a microwell, cover with the acetate sheet and incubate for 2.5 hours at room temperature (18-25°C) on an orbital microwell plate shaker with agitation (at 250 rpm).
5. Empty the contents of the microwells and wash 4 times with Wash Buffer. Washing may be performed either using microwell plate washing equipment or manually (fill the wells with Wash Buffer with a pipette or squeeze bottle, wait three minutes, empty and remove droplets by tapping the plate 4-5 times face down against absorbing material).
6. Add 100 µL of Detection Antibody to each microwell, cover with the acetate sheet and incubate the wells for 1 hour at room temperature (18-25°C) on an orbital microwell plate shaker with agitation (at 250 rpm).
7. Wash the wells by repeating Step 5.
8. Add 100 µL of Substrate to each microwell immediately after the wash step, cover the wells with the acetate sheet and incubate for 10 minutes at room temperature (18°-25°C). A blue color will develop.
9. Stop the enzymatic reaction by adding 50 µL of 0.5M H₂SO₄ to each microwell. Add the acid with the same speed and order as you added the substrate. Tap the sides of the microwell frame to ensure even distribution of the H₂SO₄. The solution color will turn yellow. Read the absorbances on a microwell plate reader at a wavelength of 450 nm within 10 minutes.


QUALITY CONTROL

The Positive Control plasma should be included each time the assay is performed. The Positive Control is diluted either 1:2 or 1:4 just prior to assaying.

RESULTS

Construct a standard curve by plotting the mean absorbance value for each standard versus its corresponding concentration of ADAMTS13. A standard curve should be generated each time the assay is performed. The following standard curve is for demonstration purposes only.

Representative Standard Curve

CALCULATIONS

The level of ADAMTS13 IgG autoantibody in the plasma sample is interpolated directly from the standard curve. If the value is greater than the highest standard, dilute the sample 1:40 or higher in Assay Buffer and repeat the assay. To obtain the level in the higher diluted sample, multiply the results by D/20 where D is the dilution factor used.

$$[\text{Autoantibody}]_{\text{Plasma Sample}} = [\text{Autoantibody}]_{\text{Diluted Test Sample}} \times D/20$$

LIMITATIONS OF THE PROCEDURE

Platelet contamination in plasma samples will interfere with the assay results. Plasma samples must be free of platelets in order to have a valid result. Exercise great care in minimizing disruption of the platelet pellet while recovering the platelet poor plasma. Samples should not be frozen and thawed more than two times.

Samples should not be collected with EDTA as the anticoagulant. Icteric, lipemic and hemolyzed samples may interfere with the assay.

EXPECTED VALUES

In a study of 22 normal plasmas (n=22), the mean level of ADAMTS13 IgG autoantibodies was found to be 4.8 AU/mL with a cut-off for positivity of 9.6 AU/mL.⁸ The cut-off for ADAMTS13 autoantibody positivity was defined as the mean + 2 SD. Each laboratory should establish its own normal range using the local population.

PERFORMANCE CHARACTERISTICS

Precision

The intra-assay coefficient of variation (CV) for this ELISA has been found to be 6.7%.

Specificity

The detection antibody used is highly specific for human IgG.

TRACEABILITY OF CALIBRATORS AND CONTROL MATERIAL

Information on traceability of calibrators and control material is available upon request.

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8. Unpublizierte Daten von American Diagnostica Inc.