american diagnostica inc.

500 West Avenue, Stamford, CT 06911-0215  
Tel. (203) 602-7777 • Fax. (203) 602-2221

# SPECTROLYSE<sup>®</sup> Blood Collection Tubes Blutabnahmeröhrchen Tubos de recolección de sangre

**REF** 452SCD

**EC** **REP**

american diagnostica GmbH  
Kaplaneigasse 35, D-64319 Pfungstadt, Germany  
Tel. + 49 61 57 / 99 08 99 • Fax + 49 61 57 / 99 08 08

ENGLISH

## INTENDED USE

SPECTROLYSE<sup>®</sup> Blood Collection Tubes are intended for the collection of whole blood for further preparation and assaying in special tests in which anti-coagulants and protease inhibitors are required to preserve the integrity of the sample. These tubes are for *in vitro* diagnostic use.

## DESCRIPTION

SPECTROLYSE blood collections tubes are non-sterile, evacuated glass tubes that contain lyophilized anti-coagulants and inhibitors specifically formulated to rapidly quench protease activity and minimize *in vitro* artifacts when assaying hemostatic parameters. The anti-coagulants and inhibitors are carefully prepared for rapid dissolution and extended shelf life. The tubes are evacuated and stoppered under controlled conditions to ensure that each tube will draw the proper volume of blood. Tubes are formulated to yield the following inhibitor concentrations in whole blood:

REF	Formulation
452SCD1	25 µM PPACK, 200 KIU/mL Aprotinin, 4.5 mM EDTA
452SCD2	25 µM PPACK, 11 mM Sodium Citrate
452SCD3	14.4 µg/mL Heparin, 15.2 µM εACA
452SCD4	75 µM PPACK

These particular anticoagulants and inhibitors were selected based upon their respective properties:

PPACK (Phe-Pro-Arg-chloromethyl ketone)	thrombin inhibition
Aprotinin	plasmin inhibition
EDTA	calcium chelation
Sodium Citrate	calcium chelation
Heparin	thrombin inhibition via AT-III
εACA (ε-Aminocaproic Acid)	fibrinolysis inhibition

## APPLICATIONS

Blood samples collected using SPECTROLYSE Blood Collection Tubes are suitable for assaying hemostatic parameters using the following formats: ELISA, amidolytic and radioimmunoassay. It is not recommended to use SPECTROLYSE tubes for clot-based assays unless the assay protocol is specifically designed to accommodate the use of these tubes.

Hemostatic Parameter	452SCD1	452SCD2	452SCD3	452SCD4
-antiplasmin	No	Yes	ND	Yes
D-Dimer	ND	Yes	ND	ND
Fibrinogen (See Reference 1)	Yes	Yes	No	Yes
Fibrinopeptide A (FPA)	Yes	No	ND	No
Heparin (via Anti-Xa)	Yes	Yes	ND	ND
Plasminogen	No	Yes	ND	ND
Plasmin- -antiplasmin Complex (PAP)	Yes	No	ND	No
PAI-1	No	Yes	ND	ND
Platelet Aggregation	No	No	ND	Yes
Prothrombin Fragment 1+2 (· <sub>2</sub> )	Yes	Yes	ND	ND
Thrombin-Antithrombin (TAT)	Yes	Yes	No	ND
TAFI Activity	No	No	Yes	No
tPA	No	Yes	ND	ND

ND – Not Determined

## USAGE/ EQUIPMENT REQUIRED BUT NOT PROVIDED

Sterile blood collection/infusion set with 12 inch tubing and multi-luer adapter (e.g. BD Vacutainer® REF 364919)  
Sterile syringe.

## USAGE/SPECIAL PRECAUTIONS

SPECTROLYSE Blood Collection tubes are *not sterile* and should *not be used for blood collection via direct phlebotomy*. Blood may be collected using a sterile collection/infusion set with a catheter of at least five inches, and equipped with a multi-sample luer adapter to eliminate the possibility of back-flush from the non-sterile tube to the patient. Alternately, blood specimens may be collected using a sterile syringe and then drawn into SPECTROLYSE tubes.

Follow the policies and procedures of your laboratory for plasma preparation and handling of potentially infectious specimens.

Blood-to-additive variation can cause erroneous results if the tubes are over filled or under filled.

Storage of glass tubes containing blood at or below 0°C may result in tube breakage.

## PRESENTATION

SPECTROLYSE Blood Collection Tubes are available in boxes of 50 tubes in 3 different fill sizes: 3 mL, 5 mL, and 10 mL tubes. When ordering, cite the REF number for the desired tube formulation followed by ordering extension denoting the desired tube size.

Formulation	Ordering No. With Tube Size		
	3 mL Tube	5 mL Tube	10 mL Tube
PPACK, Aprotinin, EDTA	452SCD1/3	42SCD1/5	452SCD1/0
PPACK, Sodium Citrate	452SCD2/3	452SCD2/5	452SCD2/0
Heparin, εACA	NA	452SCD3/5	NA
PPACK	452SCD4/3	452SCD4/5	452SCD4/0

## STORAGE AND STABILITY

Store SPECTROLYSE tubes at 2°-8°C. Tubes may be used until the recorded expiration date.

## REFERENCES

1. Tracy, R. P., D. Z. Rubin, et al. (1997). "Thrombolytic Therapy and Proteolysis of Factor V." *J Am Coll Cardiol* 30(3): 716-24.
2. Granger, C. B., J. M. Miller, et al. (1995). "Rebound Increase in Thrombin Generation and Activity after Cessation of Intravenous Heparin in Patients with Acute Coronary Syndromes." *Circulation* 91(7): 1929-35.
3. Cannon, C. P., C. H. McCabe, et al. (1997). "TNK-tissue Plasminogen Activator in Acute Myocardial Infarction. Results of the Thrombolysis in Myocardial Infarction (TIMI) 10A Dose-ranging Trial." *Circulation* 95(2): 351-6.
4. Mohler, M. A., C. J. Refino, et al. (1986). "D-Phe-Pro-Arg-chloromethylketone: its Potential Use in Inhibiting the Formation of In Vitro Artifacts in Blood Collected During Tissue-type Plasminogen Activator Thrombolytic Therapy." *Thromb Haemost* 56(2): 160-4.
5. Sie, P., A. Camez, et al. (1991). "Locally Activated Coagulation Induced by a Peripheral Endovenous Catheter." *Thromb Res* 64(6): 775-82.
6. Bauer, K. A., S. Barzegar, et al. (1991). "Influence of Anticoagulants Used for Blood Collection on Plasma Prothrombin Fragment <sub>1-2</sub> Measurements." *Thromb Res* 63(6): 617-28.
7. Rahr, H. B., J. V. Sorensen, et al. (1994). "Markers of Coagulation and Fibrinolysis in Blood Drawn into Citrate with and without D-Phe-Pro-Arg-Chloromethylketone (PPACK)." *Thromb Res* 73(5): 279-84.
8. Kereiakes, D. J., T. Lorenz, et al. (2001). "Differential Effects of Citrate versus PPACK Anticoagulation on Measured Platelet Inhibition by abciximab, eptifibatid and tirofiban." *J Thromb Thrombolysis* 12(2): 123-7.
9. Tardiff, B. E., L. K. Jennings, et al. (2001). "Pharmacodynamics and Pharmacokinetics of Eptifibatid in Patients with Acute Coronary Syndromes: Prospective Analysis from PURSUIT." *Circulation* 104(4): 399-405.

Vacutainer is a registered trademark of Becton, Dickinson and Company, USA

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#### LITERATURHINWEISE

- Tracy, R. P., D. Z. Rubin, et al. (1997). "Thrombolytic Therapy and Proteolysis of Factor V." *J Am Coll Cardiol* **30**(3): 716-24.
- Granger, C. B., J. M. Miller, et al. (1995). "Rebound Increase in Thrombin Generation and Activity after Cessation of Intravenous Heparin in Patients with Acute Coronary Syndromes." *Circulation* **91**(7): 1929-35.
- Cannon, C. P., C. H. McCabe, et al. (1997). "TNK-tissue Plasminogen Activator in Acute Myocardial Infarction. Results of the Thrombolysis in Myocardial Infarction (TIMI) 10A Dose-ranging Trial." *Circulation* **95**(2): 351-6.
- Mohler, M. A., C. J. Refino, et al. (1986). "D-Phe-Pro-Arg-chloromethylketone: its Potential Use in Inhibiting the Formation of In Vitro Artifacts in Blood Collected During Tissue-type Plasminogen Activator Thrombolytic Therapy." *Thromb Haemost* **56**(2): 160-4.
- Sie, P., A. Camez, et al. (1991). "Locally Activated Coagulation Induced by a Peripheral Endovenous Catheter." *Thromb Res* **64**(6): 775-82.
- Bauer, K. A., S. Barzegar, et al. (1991). "Influence of Anticoagulants Used for Blood Collection on Plasma Prothrombin Fragment <sub>2</sub> Measurements." *Thromb Res* **63**(6): 617-28.
- Rahr, H. B., J. V. Sorensen, et al. (1994). "Markers of Coagulation and Fibrinolysis in Blood Drawn into Citrate with and without D-Phe-Pro-Arg-Chloromethylketone (PPACK)." *Thromb Res* **73**(5): 279-84.
- Kereiakes, D. J., T. Lorenz, et al. (2001). "Differential Effects of Citrate versus PPACK Anticoagulation on Measured Platelet Inhibition by abciximab, eptifibatid and tirofiban." *J Thromb Thrombolysis* **12**(2): 123-7.
- Tardiff, B. E., L. K. Jennings, et al. (2001). "Pharmacodynamics and Pharmacokinetics of Eptifibatid in Patients with Acute Coronary Syndromes: Prospective Analysis from PURSUIT." *Circulation* **104**(4): 399-405.

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1. Tracy, R. P., D. Z. Rubin, et al. (1997). "Thrombolytic Therapy and Proteolysis of Factor V." *J Am Coll Cardiol* 30(3): 716-24.
2. Granger, C. B., J. M. Miller, et al. (1995). "Rebound Increase in Thrombin Generation and Activity after Cessation of Intravenous Heparin in Patients with Acute Coronary Syndromes." *Circulation* 91(7): 1929-35.
3. Cannon, C. P., C. H. McCabe, et al. (1997). "TNK-tissue Plasminogen Activator in Acute Myocardial Infarction. Results of the Thrombolysis in Myocardial Infarction (TIMI) 10A Dose-ranging Trial." *Circulation* 95(2): 351-6.
4. Mohler, M. A., C. J. Refino, et al. (1986). "D-Phe-Pro-Arg-chloromethylketone: its Potential Use in Inhibiting the Formation of In Vitro Artifacts in Blood Collected During Tissue-type Plasminogen Activator Thrombolytic Therapy." *Thromb Haemost* 56(2): 160-4. (Aprotinin Will Interfere with Plasminogen and □<sub>2</sub> antiplasmin Assays).
5. Sie, P., A. Camez, et al. (1991). "Locally Activated Coagulation Induced by a Peripheral Endovenous Catheter." *Thromb Res* 64(6): 775-82.
6. Bauer, K. A., S. Barzegar, et al. (1991). "Influence of Anticoagulants Used for Blood Collection on Plasma Prothrombin Fragment +<sub>2</sub> Measurements." *Thromb Res* 63(6): 617-28.
7. Rahr, H. B., J. V. Sorensen, et al. (1994). "Markers of Coagulation and Fibrinolysis in Blood Drawn into Citrate with and without D-Phe-Pro-Arg-Chloromethylketone (PPACK)." *Thromb Res* 73(5): 279-84.
8. Kereiakes, D. J., T. Lorenz, et al. (2001). "Differential Effects of Citrate versus PPACK Anticoagulation on Measured Platelet Inhibition by abciximab, eptifibatide and tirofiban." *J Thromb Thrombolysis* 12(2): 123-7.
9. Tardiff, B. E., L. K. Jennings, et al. (2001). "Pharmacodynamics and Pharmacokinetics of Eptifibatide in Patients with Acute Coronary Syndromes: Prospective Analysis from PURSUIT." *Circulation* 104(4): 399-405.

Vacutainer es una marca registrada de Becton, Dickinson and Company, USA